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38% Polyamid-Klettvelours, 28% Baumwolle, 28% Polyurethan, 6% Polyester

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GB

Wrist Supports give relief in case of injured tissue structures (joint, ligaments, capsules, cartilages). The BORT ManuCarpal®-SOFT Brace holds the fingers, thumb and the wrist in a resting position. The inserted aluminium splint enables individual adjustment in thumb and wrist area. Additionally the splint offers increased protection against wrong movements.

**Instructions for use:**

- Put brace onto the volar side and decide on desired thumb and hand position. Take off brace and bend it manually into the defined shape.
- Put on brace again and first fix it with the Velcro strap over metacarpal joint, then close the straps around the wrist and lower arm.

**Indications:**

For overnight resting, carpal tunnel syndrome. Resting in the case of tendosynovitis, postoperative.

**Contraindications:**

Lymphatic drainage disorders and indistinct soft tissue swellings remote from the applied aid Hypesthesia and circulatory disorders of the arms. Skin disorders in the treated part of the body.

**Composition:**

38% Polyamide velour, 28% cotton, 28% polyurethane, 6% polyester

**Important instructions:**

No reutilisation – this medical device is meant for the care of one patient only. In case of doctor's order to wear brace during night, avoid hindering blood circulation. In case of numbness loosen brace and remove it if necessary. If complaints continue, please contact the doctor.

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